

The Freedom[®] Knee stemmed tibial components represent a true advancement in size and fit compatibility. The component sizing chart showcases the versatility of the system and the flexibility it provides surgeons in selecting the most appropriate option for each patient.

		Femoral Component Size									
		A	B	C	D	E	F	G	H		
Stemmed Tibial Base Size	1									1/2	Tibial Liner Size / Match
	2									3/4	
	3									5/6	
	4									7/8	
	5										
	6										
	7										
	8										

		Stemmed Tibial Base Size							
		1	2	3	4	5	6	7	8
Tibial Augment Size	8								
	7								
	6								
	5								
	4								
	3								
	2								
	1								
	0								



NUMBER OF AUGMENTS	SCREW
1	5mm
2	10mm
3	15mm

AUGMENT NOTES

- Augments in orange match the Stemmed Tibial Base size
- Augments can be stacked up to 3x high
- Smaller size Augments in gray can be used to create a conical profile

ADDITIONAL NOTES

- All Offset Junctions can be used with all Stemmed Tibial Baseplates
- All Stem Extensions can be used with all Stemmed Tibial Baseplates
- Stem Extensions can be used with or without an Offset Junction

STEMMED TIBIAL SIZING

COMPONENT	AP	ML	KEEL HEIGHT
Stemmed Tibial Tray 1	40	59	25
Stemmed Tibial Tray 2	40	62	25
Stemmed Tibial Tray 3	42	66	25
Stemmed Tibial Tray 4	46	66	25
Stemmed Tibial Tray 5	48	71	25
Stemmed Tibial Tray 6	50	72	25
Stemmed Tibial Tray 7	52	76	25
Stemmed Tibial Tray 8	54	78	25

COMPONENT	AP	ML	HEIGHT
Tibial Augment Size 1	38	28	5
Tibial Augment Size 2	38	29	5
Tibial Augment Size 3	39	31	5
Tibial Augment Size 4	42	31	5
Tibial Augment Size 5	45	33	5
Tibial Augment Size 6	47	34	5
Tibial Augment Size 7	49	36	5
Tibial Augment Size 8	51	37	5

COMPONENT	THREAD	TOTAL LENGTH
Augment Screw 5mm	M4 x 0.7	7
Augment Screw 10mm	M4 x 0.7	12
Augment Screw 15mm	M4 x 0.7	17

COMPONENT	OFFSET	JUNCTION HEIGHT
Offset Junction 4mm	4	29
Offset Junction 6mm	6	29

COMPONENT	STEM DIAMETER	STEM LENGTH	SLOT
Stem Extension Ø7.5x75mm	7.5	75	–
Stem Extension Ø7.5x100mm	7.5	100	–
Stem Extension Ø7.5x150mm	7.5	150	–
Stem Extension Ø9x40mm	9	40	–
Stem Extension Ø9x75mm	9	75	–
Stem Extension Ø9x100mm	9	100	–
Stem Extension Ø9x150mm	9	150	–
Stem Extension Ø10.5x75mm	10.5	75	–
Stem Extension Ø10.5x100mm	10.5	100	✓
Stem Extension Ø10.5x150mm	10.5	150	✓
Stem Extension Ø12x40mm	12	40	–
Stem Extension Ø12x75mm	12	75	–
Stem Extension Ø12x100mm	12	100	✓
Stem Extension Ø12x150mm	12	150	✓
Stem Extension Ø13.5x75mm	13.5	75	–
Stem Extension Ø13.5x100mm	13.5	100	✓
Stem Extension Ø13.5x150mm	13.5	150	✓
Stem Extension Ø15x40mm	15	40	–
Stem Extension Ø15x75mm	15	75	–
Stem Extension Ø15x100mm	15	100	✓
Stem Extension Ø15x150mm	15	150	✓
Stem Extension Ø16.5x75mm	16.5	75	–
Stem Extension Ø16.5x100mm	16.5	100	✓
Stem Extension Ø16.5x150mm	16.5	150	✓

Freedom Total Knee® System (CR and PS)

Rx only

Maxx Orthopedics, Inc.
Plymouth Meeting, PA
19462 USA

CE 0297



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee® System.

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Freedom Knee® is manufactured by Maxx Orthopedics, Inc. FREEDOM and FREEDOM KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

© 2013 Maxx Orthopedics. All rights reserved. Updated March 2013.

Maxx Orthopedics Brochure No. MXO-MP00019-R01