

Updated 6/12

The Freedom® Knee components represent a true advancement in size and fit compatibility. The component sizing chart showcases the versatility of the system and the flexibility it provides surgeons in selecting the most appropriate option for each patient.

[illegible]

NOTE

- The Freedom Total Knee® System's Size "D" femoral component acts as a transitional size between the smaller "C" and larger "E" femoral sizes. Whereas the Freedom Knee's components typically increase in the M/L dimension in 4mm increments, the "D" is only 2mm larger than the "C" and 2mm smaller than the "E" in the M/L dimension. This transitional size allows for more flexibility in sizing up or down between tibial and femoral components and provides a greater range of fit options in pairing these components.
- The above chart represents the core product sizes available for the Freedom® Knee. Additional sizes are available upon request.

COMPONENT	AP	ML
Femoral A (Universal)	51	54
Femoral B (Universal)	54	58
Femoral C (Universal)	58	62
Femoral D (Universal)	60	64
Femoral E (Universal)	62	66
Femoral F (Universal)	66	70
Femoral G (Universal)	70	74
Femoral H (Universal)	74	78

COMPONENT	AP	ML
Tibial Base 1 (Metalback)	40	59
Tibial Base 2 (Metalback)	40	62
Tibial Base 3 (Metalback)	42	66
Tibial Base 4 (Metalback)	46	66
Tibial Base 5 (Metalback)	48	71
Tibial Base 6 (Metalback)	50	72
Tibial Base 7 (Metalback)	52	76
Tibial Base 8 (Metalback)	54	78

TIBIAL TRIAL LINER	
Purple	A1-2
Yellow	B1-2
Rust	C1-2
Red	C3-4
Gray	D1-2
Green	D3-4
Light blue	E3-4
Blue	E5-6
Magenta	F3-4
Brown	F5-6
Tan	F7-8
Bright Yellow	G5-6
Orange	G7-8
Black	H7-8

PATELLA TRIAL IMPLANT	
Yellow	25
Green	28
Light blue	31
Brown	34
Red	37
Black	40

Freedom Total Knee[®] System (CR and PS)



Maxx Orthopedics, Inc.
Plymouth Meeting, PA
19462 USA



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee[®] System.

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Freedom Knee[®] is manufactured by Maxx Orthopedics, Inc. FREEDOM and FREEDOM KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

© 2012 Maxx Orthopedics. All rights reserved. Updated March 2012.

Maxx Orthopedics Brochure No. MXO-MP00014-R02