

Updated 12/18

The Freedom[®] Knee components represent a true advancement in size and fit compatibility. The component sizing chart showcases the versatility of the system and the flexibility it provides surgeons in selecting the most appropriate option for each patient.

| | | Femoral Component Size | | | | | | | | | |
|----------------------------------|----|------------------------|---|---|---|---|---|---|---|-----|---------------------------|
| | | A | B | C | D | E | F | G | H | | |
| Tibial Component Size (all-poly) | A1 | | | | | | | | | | |
| | B2 | | | | | | | | | | |
| | C2 | | | | | | | | | | |
| | C3 | | | | | | | | | | |
| | D2 | | | | | | | | | | |
| | D3 | | | | | | | | | | |
| | D4 | | | | | | | | | | |
| | E4 | | | | | | | | | | |
| | E5 | | | | | | | | | | |
| | F5 | | | | | | | | | | |
| | G7 | | | | | | | | | | |
| Tibial Tray Size | 1 | | | | | | | | | 1/2 | Tibial Liner Size / Match |
| | 2 | | | | | | | | | 3/4 | |
| | 3 | | | | | | | | | 5/6 | |
| | 4 | | | | | | | | | 7/8 | |
| | 5 | | | | | | | | | | |
| | 6 | | | | | | | | | | |
| | 7 | | | | | | | | | | |
| | 8 | | | | | | | | | | |

NOTE

- The Freedom Total Knee[®] System's Size "D" femoral component acts as a transitional size between the smaller "C" and larger "E" femoral sizes. Whereas the Freedom Knee's components typically increase in the M/L dimension in 4mm increments, the "D" is only 2mm larger than the "C" and 2mm smaller than the "E" in the M/L dimension. This transitional size allows for more flexibility in sizing up or down between tibial and femoral components and provides a greater range of fit options in pairing these components.
- The above chart represents the core product sizes available for the Freedom[®] Knee. Additional sizes are available upon request.

Note: All sizes are in millimeters.

| COMPONENT | AP | ML |
|-----------------------|----|----|
| Femoral A (Universal) | 51 | 54 |
| Femoral B (Universal) | 54 | 58 |
| Femoral C (Universal) | 58 | 62 |
| Femoral D (Universal) | 60 | 64 |
| Femoral E (Universal) | 62 | 66 |
| Femoral F (Universal) | 66 | 70 |
| Femoral G (Universal) | 70 | 74 |
| Femoral H (Universal) | 74 | 78 |

| COMPONENT | AP | ML |
|---------------------------|----|----|
| Tibial Base 1 (Metalback) | 40 | 59 |
| Tibial Base 2 (Metalback) | 40 | 62 |
| Tibial Base 3 (Metalback) | 42 | 66 |
| Tibial Base 4 (Metalback) | 46 | 66 |
| Tibial Base 5 (Metalback) | 48 | 71 |
| Tibial Base 6 (Metalback) | 50 | 72 |
| Tibial Base 7 (Metalback) | 52 | 76 |
| Tibial Base 8 (Metalback) | 54 | 78 |

| TIBIAL TRIAL LINER | |
|--------------------|------|
| Red | A1-2 |
| Orange | B1-2 |
| Yellow | C1-2 |
| Natural | C3-4 |
| Neon Green | D1-2 |
| Dark Green | D3-4 |
| Magenta | E3-4 |
| Violet | E5-6 |
| Aqua | F3-4 |
| Light Blue | F5-6 |
| Dark Blue | F7-8 |
| Rust | G5-6 |
| Brown | G7-8 |
| Black | H7-8 |

| PATELLA TRIAL IMPLANT | |
|-----------------------|----|
| Yellow | 25 |
| Dark Green | 28 |
| Light blue | 31 |
| Brown | 34 |
| Red | 37 |
| Black | 40 |

| TIBIAL LINER HEIGHTS |
|----------------------|
| 9 |
| 10 |
| 11 |
| 12 |
| 14 |
| 17 |

Freedom Total Knee[®] System (PCK and Stemmed Tibia)



Maxx Orthopedics, Inc.
2460 General Armistead Ave, Ste 100
Norristown, PA 19403 USA



RMS UK, Ltd. 28 Trinity Rd
Nailsea, Somerset BS 484NU UK



Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee[®] System.

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Freedom Knee[®] is manufactured by Maxx Orthopedics, Inc. FREEDOM and FREEDOM KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

© 2019 Maxx Orthopedics. All rights reserved. Updated March 2019

Maxx Orthopedics Brochure No. MXO-MP00014-R04